



Item

AX-2006 Duplicating Machine

Article Number

RW01001

Certification

CE Approved ISO9001:2000 Approved

Purpose

This unit is used for melting and mixing duplicating gel in dental laboratories.

Description

This AX-2006 Duplicating Machine's tank is made of stainless steel, which is corrosion resistant, and does not affect the duplicating gel. It features high and low heating modes, optional heating power according to the amount of duplicating gel. Easy operations thanks to automatic programs and touch controls.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 500W

Rotational Speed: 40rpm

Max. Amount: 3~5Kg

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd.China
中國天津。愛鑫醫療設備有限公司

Dimensions: 52cm x 37cm x 29cm (H W D)

Shipping Dimensions: 60cm x 45cm x 35cm (H W D)

Weight: 21kg

Parts Included

Machine × 1

Manual × 1

Power Cord × 1

Fuse × 2

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com