



Item

AX-4-10 Burnout Furnace (Manual)

Article Number

MW01002

Certification

ISO9001:2000 Approved

Purpose

This unit is used for the wax elimination of casting rings and preheating the crucible in dental laboratory.

Description

The AX-4-10 Burnout Furnace has all the basic functions. The refractory brick furnace wall is excellent for reducing mechanical deterioration and chemical corrosion. It features low maintenance and low price.

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 4000W

Temperature Rating: 1000°C

Chamber Materials: Fire Brick

Capacity: Large Casting Rings x 7 or Large Casting Rings x 4 and Small Casting Rings x 4

Chamber Size: 12 × 30 × 20cm (H W D)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd.China
中國天津。愛鑫醫療設備有限公司

Dimensions: 71cm x 75cm x 46cm (H W D)

Shipping Dimensions: 88cm x 83cm x 54cm (H W D)

Weight: 120kg

Parts Included

Machine×1

Manual×1

Thermocouple×1

Thermocouple Wire×1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com