



**Item**

AX-B5 Fine Blasting Unit

**Article Number**

PF01002

**Certification**

CE Approved ISO9001:2000 Approved

**Purpose**

This unit is used for the gentle finishing of metal frameworks, and for gross-removal of investment material and oxidation layer from clasps, crown and bridge, and rings etc.

**Description**

The AX-B5 Fine Blasting Units is the updated version of the AX-B3 Sandblasting Unit. The hard steel sandblasting bottles ensure the safety of the users. It also features built-in dust drawer for easy cleaning, and redesigned lid for dust free operations.

**Technical Data**

Power Supply:220V/110V 50/60Hz

Power Consumption:<8W

Working Pressure:4Kg/cm

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy@hotmai.com

Handpiece: Twin Pencil

Nozzle: Replaceable

Nozzle Materials: Carbide/Ceramic

Nozzle Diameter: 1.2mm/0.8mm

Particles Size: 120µm/50µm

Optional Accessories: Dust Extractor (Reference AX-SUPER800, AX-MX800)

Dimensions: 31cm×42cm×42cm (H W D)

Package Size: 38cm×49cm×51cm (H W D)

Weight: 29kg

---

### **Parts Included**

Manual × 1

Nozzle × 3

Lighting Tube × 1

Fuse × 2

Air Pipe × 1

Dust Pipe × 1

Aluminum oxide × 1

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com