



#### Item

AX-J1 Dental Laboratory Cutting and Polishing Lathe

#### Article Number

RW03001

#### Certification

ISO9001:2000 Approved

#### Purpose

The unit is used for polishing castings in dental laboratories, e.g., alloys (Co-Cr, Ni-Cr) stainless steel, flexible dentures, bridges and porcelain teeth.

#### Description

AX-J1 Dental Laboratory Cutting and Polishing Lathe features low noise, high rotational speed, easy operations, and it is equipped with a safety dust collector hood which can be connected to a dust collector for dust-free operations. It also features variable speed motor.

#### Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

### Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 250W

Rotational speed: 1~3000rpm (Stepless speeds)

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 25cm × 50cm × 20cm (H W D)

Shipping Dimensions: 32cm × 68cm × 29cm (H W D)

Weight: 13kg

---

### Parts Included

Machine

Manual × 1

Power Cord × 1

Suction Hood × 1

Corundum Disc × 1

Polishing Canvas Wheel × 1

Polishing Wheel × 1

### Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com