



**Item**

AX-J5 Dental Laboratory Polishing Compact Unit (Vacuum)

**Article Number**

RW03005

**Certification**

ISO9001:2000 Approved

**Purpose**

The unit is used for polishing castings in dental laboratories, e.g., alloys (Co-Cr, Ni-Cr) stainless steel, flexible dentures, bridges and porcelain teeth.

**Description**

AX-J5 Dental Laboratory Polishing Compact Unit features low noise, high rotational speed, easy operations, and it is equipped with a suction box with light which can be connected to a dust collector for dust-free operations. The lathe, lighting system, safety eye shield and suction system make it convenient to use.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 370W

Rotational speed: 3000rpm

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com



Aixin Medical Equipment Co.,Ltd.China  
中國天津。愛鑫醫療設備有限公司

---

Max. Pressure:26kPa

Airflow Rate:172m<sup>3</sup>/h

Max. Vacuum:-16kPa

Max. Noise:<60dB

Dimensions: 46cm × 60cm × 50cm (H W D)

Shipping Dimensions: 53cm × 68 cm × 50cm (H W D)

Weight: 48kg

---

### **Parts Included**

Machine

Manual × 1

Power Cord × 1

Fuse × 2

Polishing Canvas Wheel × 2

Lighting Tube × 1

Thin Spindle × 2

Filter Bag × 2

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com