



**Item**

AX-KZ Vacuum Former

**Article Number**

RW07001

**Certification**

ISO9001:2000 Approved

**Purpose**

This unit is a necessary dental equipment for forming various kinds of plastic sheets.

**Description**

AX-KZ Vacuum Former be installed a very good vacuum motor, The case of machine is cast aluminium materials. And this machine very safe and easy to use.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 1000W

Dimensions: 25cm x 25cm x 35cm (H W D)

Shipping Dimensions: 34cm x 25cm x 46cm (H W D)

Weight: 10kg

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com



Aixin Medical Equipment Co.,Ltd. China  
中國天津。愛鑫醫療設備有限公司

---

#### **Parts Included**

Machine×1

Manual×1

Plastic Film×2

Handle×2

#### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com