



**Item**

AX-MTA Wet Model Trimmer

**Article Number**

MP03001

**Certification**

ISO9001:2000 Approved

**Purpose**

The machine is widely used in restorative dentistry, it is a necessary tool for trimming plaster models.

**Description**

AX-MTA Wet Model Trimmer of ABS engineering plastics, it has the better insulation result. The angle can be adjusted to more than 30° which enables the users to achieve perfect trimming quality.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 550W

Rotational Speed: 2,800RPM

Water Pressure: 0.5~2kg

Disc Size: 25.5cm/2.55cm/0.65cm (OD/ID/T)

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd. China  
中國天津。愛鑫醫療設備有限公司

---

Disc Materials: Corundum/Diamond

Dimensions: 32cm x 20cm x 43cm (H W D)

Shipping Dimensions: 38cm x 36cm x 49cm (H W D)

Weight: 20kg

---

### **Parts Included**

Machine

Manual×1

Power Cord×1

Dust Pipe ×1

Drain Pipe×1

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com