



Item

AX-MTB Model Arch Trimmer

Article Number

MP04001

Certification

ISO9001:2000 Approved

Purpose

The Machine is a necessary equipment for trimming the sides of plaster models.

Description

AX-MTB Model Arch Trimmer equipped with high-quality motor, speed uniform, simple operation and rugged durability.

Technical Data

Power Supply:220V/110V 50/60Hz

Power Consumption:150W

Motor rotational speed:2,800rpm

Grinding head:C 625 M06

Optional Accessories:Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions:25cm×21cm×26cm (H W D)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd. China
中國天津。愛鑫醫療設備有限公司

Shipping Dimensions: 33cm×29cm×32cm (H W D)

Weight:8.5kg

Parts Included

Machine

Manual×1

Power Cord×1

Fuse×1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com