



**Item**

AX-P3 Recyclable Sandblaster

**Article Number**

PF02002

**Certification**

CE Approved ISO9001:2000 Approved

**Purpose**

This unit is used for the gentle finishing of metal frameworks, and for gross-removal of investment material and oxidation layer from clasps, crown and bridge, and rings etc.

**Description**

AX-P3 Recyclable Sandblaster update version of AX-P2 Recyclable Sandblaster. The redesigned lid ensure dust free operations.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: <20W

Working Pressure: 6Kg/cm

Nozzle: Replaceable

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd. China  
中國天津。愛鑫醫療設備有限公司

---

Nozzle Materials: Carbide

Nozzle Diameter: 3mm

Particles Size: 250µm

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 53cm × 33cm × 40cm (H W D)

Shipping Dimensions: 60cm × 48cm × 41cm (H W D)

Weight: 24kg

---

### **Parts Included**

Machine

Manual×1

Power Cord × 1

Nozzle × 2

Fuse × 2

Air Tube × 1

Dust Pipe × 1

Corundum × 1

### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com