



**Item**

AX-PT1 Pressure Polymerizer

**Article Number**

RW06001

**Certification**

ISO9001:2000 Approved

**Purpose**

This unit is a combination of dental press and polymerization unit.

**Description**

AX-PT1 Pressure Polymerizer skips the complicated working steps and makes working experience much more easier. With this machine, clamps for polymerization unit and dental press will no long be needed. This durable machine will not have problem if used with proper care, therefore the cost of manufacturing is lowered.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Power Consumption: 600W

Motor Power: 6W

Rated pressure: 2000kg~3500kg

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com



Aixin Medical Equipment Co.,Ltd. China  
中國天津。愛鑫醫療設備有限公司

---

Dimensions: 40cm x 28cm x 14cm (H W D)

Shipping Dimensions: 46cm x 34cm x 25cm (H W D)

Weight: 10kg

---

#### **Parts Included**

Machine×1

Manual×1

#### **Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com