



Item

AX-SUPER800 Dental Vacuum Dust Extractor

Article Number

PF04001

Certification

ISO9001:2000 Approved

Purpose

This unit is for dust extraction in dental labs. It can be used with sandblasters, lathes, grinders, and workstations etc.

Description

AX-SUPER800 Dental Vacuum Dust Extractor features compact size, low noise, and high suction power. It is a portable device which is easy to operate. It can work for 4 hours continuously.

Technical Data

Power Supply:220V/110V 50/60Hz

Power Consumption:800W

Max. Pressure:26kPa

Airflow Rate:172m³/h

Max. Vacuum:-16kPa

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy1@hotmail.com



Aixin Medical Equipment Co.,Ltd. China
中國天津。愛鑫醫療設備有限公司

Max. Noise: <60dB

Diameter of Dust Pipe: 4cm

Dimensions: 30cm×26cm×36cm (H W D)

Package Size: 35cm×31cm×41cm (H W D)

Weight: 16kg

Parts Included

Machine × 1

Manua l × 1

Power Cord × 1

Fuse × 2

Filter Bag × 1

Dust Pipe × 1

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com