



Item

AX-YT1 Dental Workstation (Triple)

Article Number

MP71004

Certification

CE Approved ISO9001:2000 Approved

Purpose

This unit is used in dental laboratories for removal of dust resulting from machining.

Description

AX-YT1 Combined Dental Workstation is suitable for dental labs with large open space. Each table is designed for three dental technicians to work in a group. The gas system on the left of the table helps technicians to work conveniently with waxes. The efficient, low noise dust collection system on the right helps technicians to grind, polish etc. This workstation also features: stainless steel plate, elbow rests, retrievable air gun, strong dust collector, gas tap and tube, large size drawers, etc. The swing arm light enables users to adjust the direction of the light for brighter working environment.

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Technical Data

Power Supply: 220V/110V 50/60Hz

Power Consumption: 800W

Air Volume: 172m³/h

Max. Pressure: -16kPa

Noise Level: <68dB

Dimensions: 85cm x 60cm x 120cm (H W D)

Shipping Dimensions: 102cm x 71cm x 127cm (H W D)

Weight: 120kg

Parts Included

Fluorescent Lighting System×3

Air Nozzle×3

Arm Rest×6

Backup Power Socket×3

Dust Extraction System×3

Tools Drawer×12

Storage Drawer×21

Pan×9

Acrylic Resin Tabletop×3 (Optional)

Hettich Drawer Slide (Germany)×66 (Optional)

Gas System×3

Technician Stool×3 (Optional, please refer to C450)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy@hotmai.com