



**Item**

EM-LT3 Simplificationl Workstation (Single)

**Article Number**

MP07203

**Certification**

CE Approved ISO9001:2000 Approved

**Purpose**

This unit is used in dental laboratories for removal of dust resulting from machining.

**Description**

The EM-LT3 Workstation features storage drawers, arm rests, and piping for dust extraction system. Dust collector is optional.

**Technical Data**

Power Supply: 220V/110V 50/60Hz

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 85cm x 76cm x 52cm (H W D)

Shipping Dimensions: 102cm x86cm x 71cm (H W D)

Weight: 50kg

**Overseas Marketing:**

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

### Parts Included

Fluorescent Lighting System × 0

Air Nozzle × 0

Arm Rest × 2

Backup Power Socket × 0

Dust Extraction System × 0

Tools Drawer × 3

Storage Drawer × 0

Pan × 0

Acrylic Resin Tabletop×1 (Optional)

Hettich Drawer Slide (Germany)×12 (Optional)

Gas System×1 (Optional)

Technician Stool×1 (Optional, please refer to C450)

### Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxy@hotmai.com