



Item

EM-LT4 Simplificationl Workstation (Double)

Article Number

MP07203

Certification

CE Approved ISO9001:2000 Approved

Purpose

This unit is used in dental laboratories for removal of dust resulting from machining.

Description

The EM-LT4 Workstation features storage drawers, arm rests, and piping for dust extraction system. Dust collector is optional.

Technical Data

Power Supply: 220V/110V 50/60Hz

Optional Accessories: Dust Extractor (Reference AX-SUPER800,AX-MX800)

Dimensions: 85cm x 152cm x 52cm (H W D)

Shipping Dimensions: 102cm x 162cm x 71cm (H W D)

Weight: 80kg

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384

Tel: +86 (0)22 23859319/23859320 (English Service)

Fax: +86 (0)22 23859321

E-mail: richard@tjax.com.cn, tjaxyl@hotmail.com

Parts Included

Fluorescent Lighting System × 0
Air Nozzle × 0
Arm Rest × 2
Backup Power Socket × 0
Dust Extraction System × 0
Tools Drawer × 3
Storage Drawer × 0
Pan × 0
Acrylic Resin Tabletop×1 (Optional)
Hettich Drawer Slide (Germany)×12 (Optional)
Gas System×1 (Optional)
Technician Stool×1 (Optional, please refer to C450)

Overseas Marketing:

Add: 6B-101, No.15 Rongyuan Road, Huayuan New-Tech Industrial Park, Tianjin, China, 300384
Tel: +86 (0)22 23859319/23859320 (English Service)
Fax: +86 (0)22 23859321
E-mail: richard@tjax.com.cn, tjaxy@hotmai.com